

### Monthly

- ⑥ Gently clean the battery compartment with a soft brush.
- ⑥ Change the wax guard. If you are unsure how to do this ask your audiologist.

### When required

- ⑥ Change the battery. Your hearing aid may beep to warn you that the battery is low.
- ⑥ Carry spare batteries with you when you go out.
- ⑥ Have your ears checked for wax by your Audiologist, Doctor or Ear Nurse Specialist. We suggest doing this every six months.



# Caring for your In-the-Ear Hearing Aid



*“Providing a higher quality of life  
through better hearing”*

**HEARING PROFESSIONALS  
NELSON MARLBOROUGH**

24 Nile Street • Nelson  
Phone 03 548 2323  
Fax 03 548 2324

54 Scott Street • Blenheim  
Phone 03 578 6688  
Fax 03 578 6678

Freephone 0800 23 23 23  
Email: [info@hearingprofessionals.co.nz](mailto:info@hearingprofessionals.co.nz)  
Web: [www.hearingprofessionals.co.nz](http://www.hearingprofessionals.co.nz)




# CARING FOR YOUR IN-THE-EAR HEARING AID

**To ensure your In-the-Ear hearing aid always performs at its best, follow these simple steps:**

## Every day

When the hearing aid is not in use:

- 1 Turn the hearing aid off.
- 2 Wipe the hearing aid with a soft dry cloth.



- 3 Remove the battery.



- 4 Place the hearing aid (with the battery compartment open) in your dehumidifier.
- 5 Screw the lid back on.



## In the morning

Before putting the hearing aid in your ear:

- 6 Check the sound outlet for wax and gently remove with a soft brush.



- 7 Check for debris on the hearing aid casing, switches and microphone screen.
- 8 Gently remove any debris with a soft brush.

## Weekly

Maintain your dehumidifier:

- 9 Remove the pillow from your dehumidifier.
- 10 Place the pillow in the microwave as shown.



- 11 Microwave the pillow for 30 seconds on HIGH.
- 12 Allow the pillow to cool; the crystals should be dark blue. If the crystals are still a pink or violet colour place the pillow in the microwave for a further 30 seconds.
- 13 Repeat as many times as necessary for the crystals to turn dark blue.
- 14 Make sure you allow the pillow to cool between each microwave session.
- 15 NEVER put hearing aids into the microwave.